

COMPOSITION:

Calcium hydroxide and silicon, iron, aluminium, sodium, potassium, bismuth, magnesium oxides, calcium phosphate.

INTENDED USE

Device is intended for use during dental treatment as material for filling and remineralisation of root canals.

When mixed with water it sets in 6-10 minutes.

MTA+ is especially for use in the following cases

- demineralization of lesions of the root canal,
- reverse filling of the apex of the root,
- direct pulp capping,
- pulp amputation,
- in the treatment of teeth with incomplete root development.

Product has remineralising effect on the hard tissues of tooth, securely seals the root canal.


MTA+ should not be used as a separate or final root canal filling.

INSTRUCTIONS FOR USE

Place on the mixing plate the content of *MTA+ powder* from the glass vial. Add 1-2 drops of *MTA+ liquid* using the pipette.

Mix them within 30 sec until consistency of the wet sand is achieved. If the consistency of the mixture is too dense or crumbly then add another drop of the liquid (or distilled water).

Insert the preparation into selected field by means of the amalgam carrier.



Use the preparation only when sharp symptoms are gone as the acid pH interferes with the setting process of MTA+ .
MTA+ powder is a dose for single use. When mixed to *MTA+ liquid* (or distilled water), it should be used within 10 minutes.
After that time the product is not suitable for use due to the progressive setting process and hardening of the mixture.
Use the package of *MTA+ liquid* within one year after first opening.

WARNING

Do not use the product on patients who are hypersensitive to the product's components.

CAUTION

Product may cause eye damage. It is irritating to mucous membrane, skin and eyes.

In case of contact rinse immediately with plenty of water and seek medical advice.

In case of consumption do not regurgitate or induce vomiting. Drink plenty of water and seek medical advice.

During the treatment it is recommended that rubber-dam, single-use gloves and protective eye-glasses are used.

STORAGE

Store in original packaging at a temperature below 25 °C. Protect from moisture.

Store out of reach of children.

Expiry date is placed on the main container.

HANDLING THE PRODUCT CONTAINER

Empty containers should be disposed of responsibly.

AVAILABLE CONTAINERS

- MTA+ MINI - 3 pieces of powder *MTA+ powder* (0,14 g) + 1 piece of liquid *MTA+ liquid* (1ml)

- MTA+ STANDARD - 6 pieces of powder *MTA+ powder* (0,14 g) + 1 piece of liquid *MTA+ liquid* (1ml)

- MTA+ MAXI - 10 pieces of powder *MTA+ powder* (0,14 g) + 1 piece of liquid *MTA+ liquid* (1ml)

- MTA+ PRO - 10 pieces of powder *MTA+ powder* (0,14 g) + 1 piece of liquid *MTA+ liquid* (1ml) + MATRIX MTA+



WARNING. READ ENCLOSED INSTRUCTIONS CAREFULLY



MAXIMUM ALLOWABLE TEMPERATURE UP TO 25°C



PROTECT FROM MOISTURE



DO NOT RE-USE



PPH CERKAMED Wojciech Pawłowski
37-450 Stalowa Wola, Polska
ul. Kwiatkowskiego 1
Tel/fax 0048 15 842-35-85