
Dental training kit

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USER MANUAL

General Rules

ADSM201 is a multifunctional digital microscope. It allows object observation in Full HD, with high magnification, it has multifunction menu and multi-output at the same time. It's design combines many ergonomic aspects, which provide users with better experience.

- This manual, guides how to work safely when using microscope.
- Not following of this instruction may result in permanent microscope damage.
- Storage and transport: Incorrect transportation may result in device damage. Always transport device in original packaging, to prevent any damages.
- Cleaning: in case when lens cleaning is required, use appropriate liquid and cloth intended for this purpose. Clean other elements with damp cloth.

WARRANTY: The warranty does not cover defects caused by failure to follow of this instruction. The device is covered by 1-year warranty from the date of purchase. The microscope ought to be used with original accessories kit with power supply described in instruction. Do not fix, dismount or alter device on your own. It is recommended to use device in temperature from -10 to 60 degrees of Celsius.

Safety Rules


1. Store microscope in dry and clean place. Microscope is intended for indoor use only.
2. Keep kids and third party people away when using microscope.


Basic Parameters:

Image sensor	3 Mega Pixels HD Sensor
Output	HDMI&USB
Video format	MOV
Magnification Ratio	Up to 300X(HDMI monitor 22 inch)
Photo resolution	12M
Photo format	JPEG
Focus Range	10mm do 150mm
Frame Rate	Max 30f/s under 600 Lus Brightness
Digital zoom	4X
Storage	Micro-SD card, up to 32G
PC support	Windows XP/7/8
Power source	5V DC
Stand size	14.5cm*7.9cm*16.5cm
Monitor size	3 inch
Accessories	Stand*1, Adapter*1, USB Cable*1, HDMI Cable*1
Weight	1.2kg

Parts of a Digital Microscope ADSM201

General use:

1. Connect the camera to the power supply via USB cable. Also connect spotlights to the power supply.
2. Press the "ON/OFF" button on the monitor. (Sometimes the camera will switch on automatically when it is connected to the power supply)
3. Turn on the spotlights, adjust the "height adjustment" and the "Focus controller" to get the image that you want.
4. Use the "UP/DOWN" to do digital zooming.

Menu setup

1. Press "Menu" until it shows the menu on the built-in screen.
2. Press "UP/DOWN" and "OK" choose what to set.
3. Press "Menu" again after you finish setting.

To switch mode

1. Press "Mode" shortly.
2. Icons, showing actual work mode:

 : standby mode

 : capture mode

 ● : record mode

 : image playback mode

 : video playback mode – (If the files are videos, press "OK" to play the video file)

 ● 🔑 : lock the file which is being recording

To capture pictures

Press "Mode" shortly to get into the Capture mode. Press "OK" to take pictures.

To record a video

1. Press "Mode" shortly to get into the Standby mode. Press "OK" to start recording.
2. Press "OK" again to stop recording.


Displaying on an External HDMI monitor:

- First, make sure that the monitor connected to the microscope is literally HDMI monitor. Otherwise it won't perform as good as we expect.
- Second, connect the microscope to the HDMI monitor, also, connect it to the power supply. Because there is no built-in battery in the device, so make sure the microscope is always charging when you use it.
- Third, turn on the microscope and adjust the distance. Then put the object on the base and adjust the focus controller.
- It won't show anything on its own built-in screen when it is connecting to other monitors. So it's normal that when you use it as HDMI output, the 3 inch screen keeps being dark.
- You still can take pictures and save them if you have a memory card inside the microscope. But you can't record videos in this output mode.

USB output mode:

1. Connect the microscope to the computer using USB cable.
2. It will show 2 options: PC camera & Mass storage on the screen of the microscope. Choose "PC camera".
3. Open the software, the title will show the Storage path for the file.
4. Record video: Capture → start capture → OK → stop capture

5. The video will automatically be saved on the selected path, showed by title.
- Change path for recorded file: File → Set capture file
6. Take photo:
Capture → Still → enable/
capture → Still → folder (change photo path)
capture → snap (just take photo or click „Enter“ on keyboard))
(P.S. Using this USB output maximum resolution is 1080*720P)

Disposal of used electrical and electronic appliances

The meaning of the symbol on the product, its accessory or packaging indicates that this product shall not be treated as household waste. Please, dispose of this product at your applicable collection point for the recycling of electrical & electronic equipment waste. The correct disposal of this product will help save valuable natural resources and help in preventing the potential negative impact on the environment and human health, which could be caused as a result of improper liquidation of waste.

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Producer: Shenzhen Andonstar Technology Co., Ltd.