PL

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
EN

Biomaterial for filling and rebuilding of root canals

PLEASE READ THE FOLLOWING INSTRUCTIONS BEFORE USE FOR DENTAL USE ONLY

INGREDIENTS:

Bio MTA+ powder: calcium oxide, hydroxyapatite, oxides of: silicon, iron, aluminum, sodium, potassium, bismuth, magnesium, zirconium; calcium phosphate.
Bio MTA+ liquid: Ph Eur purified water, calcium catalyst.

INDICATIONS FOR USE

The product is intended for dental treatment as a material for filling and remineralization of root canals. When mixed with the *Bio MTA+ liquid*, the working life is ca. 4 min, and the compound sets completely in ca. 2 h. The high content of calcium ions remineralises the tooth tissue. The silicon and calcium compounds support tissue regeneration after wall perforation and intracanal resorption. *Bio MTA+* contains hydroxyapatite, a natural component of bones, which integrates the compound perfectly into the bone structure.

Bio MTA+ is especially recommended in:

- root canal wall perforation
- intracanal resorption
- filling of the posterior root apex
- direct coating of the pulp
- pulp amputation
- treatment of teeth with interrupted root development

Bio MTA+ is not suitable as a standalone permanent root canal filling. Once its two components have been mixed, *Bio MTA+* exhibits extremely plastic and cohesive consistence without brittling and stays on the instruments. This facilitates the application *Bio MTA+* when compared with the traditional MTA compounds.

INSTRUCTIONS FOR USE

Cut off the *Bio MTA+ liquid* pipette tip.

Move the entire contents of the *Bio MTA+ powder* glass vial with 1 to 2 drops of the *Bio MTA+ liquid* on a mixing plate. Mix for ca. 30 s until the compound reaches a consistence of soft plasticine (modelling putty). If the compound mix is too thick or too brittle, add one more drop of the *Bio MTA+ liquid*.

CAUTION! Mix the *Bio MTA+ powder* with the *Bio MTA+ liquid* only. Do not replace the liquid with water or other fluids. Place the compound at the target site with an applicator. One container of the *Bio MTA+ powder* contains a one-off dose. Once mixed with the *Bio MTA+ liquid*, consume the mixed dose within a maximum of 4 minutes. The compound cures and hardens afterwards.

CONTRAINDICATIONS

Do not use in patients who are hypersensitive to any of the ingredients. Do not use in inflammatory conditions. Acidic pH values disrupt the setting of *Bio MTA+*. Do not mix the *Bio MTA+ powder* with water or other fluids. Use the *Bio MTA+ liquid* only.

PRECAUTIONS

The product can cause eye damage. The product irritates the oral mucosa, the skin and eyes. In case of contact, rinse with plenty of water and contact a doctor immediately. If swallowed, do not induce vomiting. Drink plenty of water. Contact a doctor immediately. Use rubber dams, disposable gloves and safety goggles when working with the product.

STORAGE

Store in a tightly closed original container below 25 °C. Protect from moisture. Keep out of the reach of children. The expiry date is on the direct packaging.

DISPOSAL OF PACKAGING

Waste containers should be disposed of or returned to the manufacturer.

PACKAGING

Glass vial containing 0.14 g of the *Bio MTA+ powder*; pipette containing 1 ml of the *Bio MTA+ liquid*.

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