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
Gel promoting complete polymerisation of dental composites

PLEASE READ THE FOLLOWING INSTRUCTIONS BEFORE USE FOR DENTAL USE ONLY

INGREDIENTS: propanetriol (glycerol), Ph Eur purified water, gel base.

INDICATIONS FOR USE

INOX Gel is intended for use during:
- polymerisation process of light-hardened composites;
- prosthetic treatments performed with the use of dental cement (e.g. placing veneers, inlays, onlays), in order to promote complete polymerisation of composites.
INOX Gel protects the surface of composites by blocking the access of oxygen during the light-curing process. In this way, an inhibition layer does not form on the composite surface, so that the polymerization process is not inhibited and can proceed properly.
The composite surface hardened with INOX Gel is smooth, harder, more scratch- and abrasion-resistant. The composite is also non-sticky, resistant to discolouration and does not wash out.
The light of a polymerisation lamp penetrates through the gel layer. The product is easily washed out with water.
NOTE: When using a composite material, strictly follow the manufacturer's recommendations.
INOX can also be used as a try-in gel for cementing of prosthetic works.

INSTRUCTIONS FOR USE

Polymerisation of light-cured composites:

1. Apply and form the last composite coat, apply INOX Gel without touching the composite coat with the applicator.
- The gel layer must tightly cover the whole surface of the dental restoration and stay in place during the polymerisation process.
2. Light cure through the gel coat.
3. Wash out INOX Gel after the polymerisation process has been finished.
4. Trim and polish the composite surface as needed.

During seating of prosthetic works with resin cement:

1. When the restoration has been seated, remove excess resin cement as needed.
2. Apply the INOX Gel on the marginal area of the restoration (veneer, inlay, onlay) without touching the composite coat with the applicator. The gel layer must tightly cover the whole surface of the dental restoration and stay in place during the polymerisation process.
3. Light cure through the gel coat.
4. Wash out INOX Gel after the polymerisation process has been finished.
5. Trim and polish the composite surface as needed.

This product is reusable, but the applicator is intended for single use only. Repeated use of the applicator can pose a risk of secondary infections.

CONTRAINDICATIONS

Do not use in patients sensitive to any of the ingredients.

PRECAUTIONS

Use a dental dam, gloves, safety goggles and protective clothing when working with the product.
This product can cause eye irritation. In case of contact with eyes, rinse with plenty of water.

STORAGE


Store in tightly closed original container below 25°C.
Keep out of the reach of children.
The expiry date is on the immediate packaging.
After opening of the packaging the expiry date does not change, providing that it is tightly sealed after each use.

DISPOSAL OF PACKAGING

Used containers should be disposed of or returned to the manufacturer.

PACKAGING

A syringe containing 2 ml of the formulation, disposable applicators.

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