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EN Material for filling and rebuilding of root canals

PLEASE READ THE FOLLOWING INSTRUCTIONS BEFORE USE

TO BE USED BY DENTIST ONLY

COMPOSITION:

MTA+ powder: kalcium oxide, oxides of: silicon, iron, aluminium, sodium, potassium, bismuth, magnesium, zirconium; calcium phosphate, MTA+ liquid: aqua Ph. Eur. purified water.

INTENDED USE:

The product is intended for dental treatment as a material for filling and remineralization of root canals. The compound cures in 6 to 10 minutes after mixing with water. The high content of calcium ions remineralizes the tooth tissue. The silicon and calcium compounds support tissue regeneration after wall perforation and intracanal resorption. MTA+ is especially recommended in:

- root canal wall perforation
- intracanal resorption
- filling of the posterior root apex
- direct coating of the pulp
- pulp amputation
- treatment of teeth with interrupted root development.

MTA+ is not suitable as a standalone permanent root canal filling.

INSTRUCTIONS FOR USE:

Cut off the MTA+ liquid pipette tip. Move the entire contents of the MTA+ powder glass vial with 1 to 2 drops of the MTA+ liquid on a mixing plate.

Mix for ca. 30 seconds until the compound reaches a consistency of moist sand. If the compound is too dry add a drop of the liquid (or deionized water). Place the compound at the target site with an applicator.

Exposure to moisture, including blood, will reduce the curing time.

One container of the MTA+ powder contains a single dose.

Once mixed with the MTA+ liquid (or deionized water) the mixed dose must be used within 10 minutes. The compound cures and hardens afterwards.

CONTRAINDICATIONS:

Do not use in patients who are hypersensitive to any of the ingredients. Do not use in inflammatory conditions. Acidic pH values disrupt setting of MTA+.

PRECAUTIONS:

Causes serious eye damage. Irritating to mucous membranes of the oral cavity and skin. In case of contact with skin, wash with plenty of water and contact a doctor immediately. If swallowed, do not induce vomiting. Drink plenty of water. Contact a doctor immediately. Use rubber dams, disposable gloves and safety goggles when working with the product.

Any serious incident involving the product should be reported to the manufacturer and the competent authority of the State in which the user or patient is domiciled.

STORAGE:

Store in a tightly closed original container below 25°C.

Protect from moisture.

Keep out of reach of children.

The expiry date is on the direct packaging.

HANDLING THE PRODUCT CONTAINER:

Empty containers should be disposed of responsibly or returned to the manufacturer.

PACKAGING:

Glass vial containing 0,14 g of the MTA+ powder; pipette containing 1 ml of the MTA+ liquid.

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