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EN
Biomaterial for filling and rebuilding of root canals
**PLEASE READ THE FOLLOWING INSTRUCTIONS BEFORE USE TO BE USED BY DENTIST ONLY**

COMPOSITION:
*BIO MTA+ powder*: calcium oxide, hydroxyapatite, oxides of: silicon, iron, aluminum, sodium, potassium, bismuth, magnesium, zirconium, calcium phosphate.
*BIO MTA+ liquid*: Ph. Eur. purified water, calcium catalyst.

INTENDED USE
The product is intended for dental treatment as a material for filling and remineralisation of root canals. After mixing with the liquid the *BIO MTA+ liquid* has a working time of approx. 4 min and the material sets fully after approximately 180 min (approx. 3 h).
The high content of calcium ions remineralises the tooth tissue. The silicon and calcium compounds support tissue regeneration after wall perforation and intracanal resorption. *BIO MTA+* contains hydroxyapatite, a natural component of bones, which integrates the compound perfectly into the bone structure.
*BIO MTA+* is especially recommended in:
- root canal wall perforation
- intracanal resorption
- filling of the posterior root apex
- direct coating of the pulp
- pulp amputation
-treatment of teeth with interrupted root development.
*BIO MTA+* is not suitable as a standalone permanent root canal filling. Once its two components have been mixed, *BIO MTA+* exhibits extremely plastic and cohesive consistency without brittleing and stays on the instruments. This facilitates the application *BIO MTA+* when compared with the traditional MTA compounds.

INSTRUCTIONS FOR USE
Cut off the *BIO MTA+ liquid* pipette tip. Move the entire contents of the *BIO MTA+ powder* glass vial with 1 to 2 droplets of the *BIO MTA+ liquid* on a mixing plate. Mix for ca. 30 s until the compound reaches a consistency of soft plastincine (modelling putty). If the compound mix is too thick or too brittle, add one more drop of the *BIO MTA+ liquid*.

**CAUTION!** Mix the *BIO MTA+ powder* with the *BIO MTA+ liquid* only. Do not replace the liquid with water or other fluids.
Place the compound at the target site with an applicator.
One container of the *BIO MTA+ powder* contains a one-off dose. Once mixed with the *BIO MTA+ liquid*, consume the mixed dose within a maximum of 4 minutes. The compound cures and hardens afterwards.

CONTRAINDICATIONS
Do not use in inflammatory conditions. Acidic pH values disrupt the setting of *BIO MTA+*. The *BIO MTA+* cannot be used as a full root canal filling, as its potential removal will be difficult.
Do not use *BIO MTA+ powder* as a stand-alone permanent root canal filling.
Do not mix the *BIO MTA+ powder* with water or other fluids. Use the *BIO MTA+ liquid* only.

PRECAUTIONS
Causes serious eye damage.
Irritating to the mucous membranes of the oral cavity and skin.
May cause respiratory irritation. Exposure to moisture, including blood, reduces the curing time. Completing filling of the root canal may reduce the resistance of the root to fracture.

Discolouration might occur. Multi-stage dental treatment may involve temporary pain.
The use of permanent dental fillings can be associated with microleakage. MTA products have a cytotoxic potential. Use a cofederant, disposable gloves, safety goggles/visors and protective clothing when handling the product. In case of contact, rinse immediately with plenty of water. If swallowed, rinse the mouth with water. Do not induce vomiting. Avoid creating and inhaling dust. In case of worrying symptoms, consult a doctor.
Any serious incident involving the product should be reported to the manufacturer and the competent authority of the State in which the user or patient is domiciled.

STORAGE
Store in a tightly closed original container below 25°C.
Do not use the material after the expiry date of *BIO MTA+ liquid*.
Keep out of reach of children.
The expiry date is on the direct packaging.

HANDLING THE PRODUCT CONTAINER
Empty containers should be disposed of responsibly or returned to the manufacturer.

PACKAGING
Glass vial containing 0,14 g of the *BIO MTA+ powder*; pipette containing 1ml of the *BIO MTA+ liquid*.

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EN
Biomaterial for filling and rebuilding of root canals
**PLEASE READ THE FOLLOWING INSTRUCTIONS BEFORE USE TO BE USED BY DENTIST ONLY**

COMPOSITION:
*BIO MTA+ powder*: calcium oxide, hydroxyapatite, oxides of: silicon, iron, aluminum, sodium, potassium, bismuth, magnesium, zirconium, calcium phosphate.
*BIO MTA+ liquid*: Ph. Eur. purified water, calcium catalyst.
INTENDED USE
The product is intended for dental treatment as a material for filling and remineralisation of root canals. After mixing with the liquid the *BIO MTA+ liquid* has a working time of approx. 4 min and the material sets fully after approximately 180 min (approx. 3 h).
The high content of calcium ions remineralises the tooth tissue. The silicon and calcium compounds support tissue regeneration after wall perforation and intracanal resorption. *BIO MTA+* contains hydroxyapatite, a natural component of bones, which integrates the compound perfectly into the bone structure.
*BIO MTA+* is especially recommended in:
- root canal wall perforation
- intracanal resorption
- filling of the posterior root apex
- direct coating of the pulp
- pulp amputation
-treatment of teeth with interrupted root development.
*BIO MTA+* is not suitable as a standalone permanent root canal filling. Once its two components have been mixed, *BIO MTA+* exhibits extremely plastic and cohesive consistency without brittleing and stays on the instruments. This facilitates the application *BIO MTA+* when compared with the traditional MTA compounds.

INSTRUCTIONS FOR USE
Cut off the *BIO MTA+ liquid* pipette tip. Move the entire contents of the *BIO MTA+ powder* glass vial with 1 to 2 droplets of the *BIO MTA+ liquid* on a mixing plate. Mix for ca. 30 s until the compound reaches a consistency of soft plastincine (modelling putty). If the compound mix is too thick or too brittle, add one more drop of the *BIO MTA+ liquid*.

**CAUTION!** Mix the *BIO MTA+ powder* with the *BIO MTA+ liquid* only. Do not replace the liquid with water or other fluids.
Place the compound at the target site with an applicator.
One container of the *BIO MTA+ powder* contains a one-off dose. Once mixed with the *BIO MTA+ liquid*, consume the mixed dose within a maximum of 4 minutes. The compound cures and hardens afterwards.

CONTRAINDICATIONS
Do not use in inflammatory conditions. Acidic pH values disrupt the setting of *BIO MTA+*. The *BIO MTA+* cannot be used as a full root canal filling, as its potential removal will be difficult.
Do not use *BIO MTA+ powder* as a stand-alone permanent root canal filling.
Do not mix the *BIO MTA+ powder* with water or other fluids. Use the *BIO MTA+ liquid* only.

PRECAUTIONS
Causes serious eye damage.
Irritating to the mucous membranes of the oral cavity and skin.
May cause respiratory irritation. Exposure to moisture, including blood, reduces the curing time. Completing filling of the root canal may reduce the resistance of the root to fracture.

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